**M100 SM 1-wire/M100 SM RS485**

**Sensor Mount Transmitter**

**Instruction Manual**

–	Regulatory approval	2
cs	Návod k obsluze	8
da	Brugsanvisning	12
de	Betriebsanleitung	16
en	Instruction manual	20
es	Manual de instrucciones	24
fi	Käyttöohje	28
fr	Instructions d'utilisation	32
hu	Használati utasítás	36
it	Istruzioni per l'uso	40
ja	取扱説明書	44
ko	사용 매뉴얼	48
nl	Gebruikershandleiding	52
pl	Instrukcja obsługi	56
pt	Manual de instruções	60
ru	Инструкция по эксплуатации	64
sv	Bruksanvisning	68
th	คู่มือคำแนะนำ	72
zh	说明书	76

## Regulatory approval

### FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

### IC Statement

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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### Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

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## Regulatory approval

For regulatory approvals please refer to pages 2 to 7.

## 1 Introduction

The M100 SM is a direct mount transmitter for METTLER TOLEDO ISM® sensors (1-wire sensors or RS485 digital output sensors) for use with analog or digital biocontrollers.

The M100 SM provides Bluetooth 4.0 communication between the ISM sensor and the METTLER TOLEDO iSense™ software or iSense Mobile app. (iSense and iSense Mobile enable ISM sensor calibration, configuration, diagnostics and data logging.)

The M100 SM stores the configured span/range and has a current output that can be aligned with the stored span/range.

Plug and Measure means configuration of a sensor upon connection to an M100 SM is unnecessary.

An independent RS485 MODBUS RTU interface is also available for the M100 SM.

For communication with a process logistics controller (PLC/DCS), the M100 SM features two standardized analog current outputs 4/20 mA or with a digital RS485 digital interface (MODBUS RTU). The M100 SM is double-galvanically isolated and provides a stable analog or digital signal.

## 2 Technical data M100 SM

Measurement parameters	1-wire: pH, amperometric DO, CO <sub>2</sub> RS485: optical DO
ISM	Full ISM functionality
Transmission frequency	2.4 GHz, Bluetooth 4.0
Range	10 m
Transmission interval	1 s
Power supply	24 V (min. 375 mA)
Operating temperature	-15 to +60 °C
Relative humidity	5 ... 95 % rH (non-condensing)
4 – 20 mA current range	3.5 to 22 mA (NAMUR)
Protection class	IP67

### 2.1 Connecting sensors

The M100 SM 1-wire has a connector to fit any METTLER TOLEDO 1-wire ISM sensor (AK9 connector).

The M100 SM RS485 has a connector that fits METTLER TOLEDO's RS485 sensors (optical oxygen).

### 2.2 VP8 connector to PLC/DCS

**Note:** for connecting the M100 SM to the PLC/DCS, use METTLER TOLEDO's standard VP8 cables with open end.

PIN/cable assignment VP 8 cable:

VP pin	VP8 cable	Function
A	black/transparent	out1 4 to 20 mA+ <sup>1</sup>
B	red	out1 4 to 20 mA–
C	gray	24 VDC+
D	blue	24 VDC–
E	white	out2 4 to 20 mA+ <sup>1</sup>
F	green	out2 4 to 20 mA–
G	pink	RS485+ (MODBUS RTU)
H	brown	RS485– (MODBUS RTU)
S	green/yellow	shield

**1** The 4/20 mA output is an active source and needs to be connected to a passive receiver e.g., passive PLCard.

## 2.3 4/20 mA output configuration

The configuration of the 4/20 mA outputs can be conducted with iSense PC software or with the iSense Mobile app (download from Google Play or Apple App Store).

## 2.4 MODBUS RTU

For support of the MODBUS RTU programming please download the MODBUS Programmers Guide from [www.mt.com/pro-M100SM-documentation](http://www.mt.com/pro-M100SM-documentation).

### 3 Pairing to iSense Mobile

The interface works in conjunction with the iSense Mobile app (download from Google Play or Apple App Store) or with iSense PC software.

[www.mt.com/isense-mobile](http://www.mt.com/isense-mobile)

After activating Bluetooth and opening iSense Mobile or iSense, please select the sensor you wish to communicate with.

### 4 LED messages

LED green on	Powered up and connected to sensor
LED blue on	Bluetooth connection established with ISM sensor
LED yellow on	At least one warning registered. See diagnostics in iSense iMonitor.
LED red on	At least one alarm has been registered. See diagnostics in iSense iMonitor.

ISM is a registered trademark of the METTLER TOLEDO Group in Switzerland, Brazil, USA, China, European Union, South Korea, Russia and Singapore.

iSense is a trademark of the METTLER TOLEDO Group.

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© METTLER TOLEDO. All rights reserved
- France**                    **Mettler-Toledo Analyse Industrielle S.A.S.**
- Germany**                **Mettler-Toledo GmbH**
- Switzerland**           **Mettler-Toledo (Schweiz) GmbH**
Mettler-Toledo GmbH, Process Analytics