EasyMax™ 102

Pre-Installation Guide

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Guía de preinstalación

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Pre-Installation Guide

Site Requirements for the Standard Equipment

METTLER TOLEDO

Pre-Installation Guide

- This document outlines the site requirements for a successful installation of the EasyMax™ system and provides information about technical data, construction materials and standard delivery.
- For more detailed information please check the Product Catalog and the Operating Instructions of the EasyMax™.

Space requirements

The space requirements take into account only the physical size of the EasyMax™ including cooling connector and touchpad (USB memory stick excluded). Additional space is needed for accessories (pH meter, flow meter, cryostat etc.)

Power requirements

Voltage: 100 – 240 VAC
 Max. voltage fluctuation: ±15%
 Frequency: 50/60 Hz
 Power consumption: max. 1000 VA

Electrical connectors


- Power socket
- RS232 interface
- Ethernet
- CAN Bus
- USB
- Connector for touchpad

Weight

including touchpad: 15 kg

Cooling media

- tap water, ethanol, ethylene glycol or a silicone oil

Purge gas

In order to prevent internal corrosion due to condensation of atmospheric moisture, the EasyMax™ has to be purged with dry gas (dry air, nitrogen, argon).

| | |
|-----------------------------------|-------------------------------------------------------------------------------------|
| Min. inlet pressure of purge gas: | 0.05 bar or 0.73 psi |
| Min. gas flow: | 3 L to 5 L/min |
| Gas connections: | Quick connect couplings for tubing with 4 mm inner diameter (included in delivery). |

Notice: Operation without purge can damage the instrument and will void the warranty.

Cooling

The EasyMax™ requires a constant flow of coolant liquid. Either tap water (salt free) or a cryostat may be used. If the water is polluted or has a high degree of hardness, filters should be installed in the inlet line.

Temperature range of the cooling liquid: -20 °C to 40 °C (inlet temperature).

| | |
|------------------------|------------------------------------------------------------------------------------|
| Min. flow rate: | 3 L/min |
| 2 Coolant connections: | in/out with hose nipple for tubing with 8 mm inner diameter (included in delivery) |

If you use a cryostat for cooling, mind the following specifications (one cryostat per EasyMax™):

| | |
|------------------------|---------|
| Min. pump capacity: | 0.5 bar |
| Cooling power at 0 °C: | 350 W |

You can also use one cryostat for cooling several EasyMax™ systems (in parallel), but then you need a pump with higher capacity (min. 0.5 bar).

The cooling power requirements at 0 °C increase with every EasyMax™ connected (1: 350 W, 2: 700 W etc.).

Notice: Operation without coolant flow triggers an emergency program (temperature of the coolant outlet $T_c > 40$ °C).

Ambient conditions

According to EN 61010-1, the following requirements must be fulfilled for safe operation of the EasyMax™ system:

| | |
|-------------------------------------|------------|
| Temperature: | 5 to 40 °C |
| Max. relative atmospheric humidity: | 80% |

Material used for construction

| | |
|-------------------------------------------------------|------------------------------------------------|
| Cover plate: | stainless steel, PFA/FEP coated |
| Housing material: | powder coated stainless steel |
| Reactor windows: | Borosilicate glass 3.3 |
| Receptacles for thermostat block: | anodized aluminum |
| Sealing rings for thermostat block: | PTFE/25% carbon |
| Holder for lab bars: | stainless steel |
| Connectors for purge gas: | stainless steel, nickel-plated brass |
| Purge gas lines: | PVC, FEP, PP, PVDF, PTFE, aluminum |
| Coolant tubings: | PVC, PVDF, anodized aluminum |
| Coolant connectors: | nickel-plated brass |
| Touchpad: | PA12, aluminum |
| Protective cover for touchpad: | Barex® resin |
| Glass Reactors: | Borosilicate glass 3.3 |
| Magnetic stirrer: | PTFE coated |
| or overhead stirrer: | borosilicate glass or Hastelloy® C22 |
| LEMO connectors for Tr sensors and overhead stirrers: | chrome-plated brass with protection cap in POM |
| USB connector: | stainless steel with protection cap in POM |
| On/Off switch: | stainless steel |

Standard equipment

The EasyMax™ will be shipped with*:

- Lab bars, diameter: 14 mm
- Touchpad with 2 protective covers (cable length: 1 m)
- Pen for touchpad
- Two 2.5 m PVC tubes for coolant, inner diameter: 8 mm
- 2 m PVC tube for purge gas, inner diameter: 4 mm (incl. reduction piece)
- Country specific power cable (3 m)
- USB memory stick for data export
- Filter for tap water
- Flow controller for cooling water
- Installation Instructions and Quick Start Guide (languages: EN, DE, FR, ES, JP, CN)
- Factory Test Report (EN)
- Installation Check List (EN)
- Operating Instructions (EN)
- CE certificate

* Glassware and other accessories according to customers order.

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©06/2008 Mettler-Toledo AG
Mettler-Toledo AutoChem Inc.
Mettler-Toledo AG, AutoChem