METTLER TOLEDO

pH ISM Qualification Kit

**Instruction manual
Bedienungsanleitung
Instructions d'utilisation**

pH ISM Qualification Kit
52 300 439

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Instruction manual

Contents

	Page
1 Safety	5
2 Product Description	7
3 Operation	8
3.1 Connecting the transmitter to the pH ISM Qualification Kit	8
4 Product Specification	10
4.1 Simulator, Pre-registered Values	10
4.2 Sensor Wear Data.....	10
4.3 Warnings and Alarms	11
5 Installation	12
5.1 Dimensions	12
6 Ordering Information	13
7 Certificates	14
8 Protocol	15

1 Safety

This instruction manual lists the main points of how to use the **pH ISM Qualification Kit** from METTLER TOLEDO for best results and in the intended manner.

This **instruction manual**, and the **safety requirements** in particular, **must be observed** by the **operator and all other persons** working with this equipment.

The instruction manual must always be stored close at hand, in a place accessible to all people working with an ISM pH sensor and the pH ISM Qualification Kit.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of product.

Intended use

The pH ISM Qualification Kit is intended **solely for validating and testing of the ISM transmitter family**, respectively, as described in this instruction manual.

The pH ISM Qualification Kit is intended **for simulating the functions of METTLER TOLEDO ISM sensors**, as described in this instruction manual.

Any other use, or use not mentioned here, that is incompatible with the technical specifications is deemed inappropriate. The operator is solely responsible for any damage arising from such use.

Other prerequisites for appropriate use include:

- observing the instructions, notes and requirements set out in this instruction manual
- correct care and maintenance of the unit, according to the instruction manual
- observing the prescribed environmental and operational conditions
- compliance with local statutes

2 Product Description

The pH ISM Qualification Kit is a set of 5 different sensor plug-in heads and allows the simulation of reading values of pH ISM sensors with predefined values and errors (not changeable by the user) to assure safe use as qualification tool.

Each sensor plug-in head corresponds to one of METTLER TOLEDO's pH ISM sensors and delivers the complete data of an ISM sensor.

The pH ISM Qualification Kit has basically 2 functions:

- Tool for the qualification of measuring systems.
- Service tool for quick checks and loop trouble-shooting.

The data indicated (primary/secondary) and transmitted are as follow:

- Main measured value
- Temperature
- Error message according to indication.

All standard process and wear data of the sensor are also available in each plug-in head and are also sent to the transmitter or the PLC as well as to the iSense software.

3 Operation

3.1 Connecting the transmitter to the pH ISM Qualification Kit

To qualify your measuring systems proceed as follows:

Please refer to the transmitter instruction manual for the exact description of the transmitter operation.

Step 1: Disconnect the ISM sensor from the measuring point.

Step 2: Connect one of the pH ISM Qualification Kit plug-in heads to the transmitter via the standard connecting cable.

Step 3: The preregistered values then are displayed on the transmitter.

Step 4: Check if the displayed values by the transmitter or the PLC correspond to the predefined fixed values of the pH ISM Qualification Kit plug-in head.

 See certification for exact values and tolerances.

Step 5: In case you are using the Error simulators, check that the appropriate errors are correctly displayed by the transmitter or the PLC.

Step 6: Report results into the attached spreadsheet "Protocol for pH ISM Qualification Kit".

Problem solving: If any of the tests described above failed, you should change the connecting cable and repeat the test. If the test still shows the same result, the transmitter or the cable is defective. If you are using a transmitter with current outputs, so ensure the outputs are correctly wired.

You must **never do a calibration** (the simulators are read-only).

You must **never increment autoclave cycles** (always set "NO" to cycle counter), this is only needed if autoclave has been set to a value higher than 0.

Different transmitters may display different values and errors.


Out of range behavior is also different depending on the transmitter in use thus **different values/readings may occur.**

Set the transmitter input filter to off (only for SIM pH TGL).

Do not use the pH ISM plug-in heads after the expiry date shown on the label of each plug-in head. Please send them back to METTLER TOLEDO for recertification.

Measuring point with
METTLER TOLEDO Electrode

Connecting pH ISM
Qualifying plug-in head

Environment operating and storing conditions:

Permissible ambient temperature (operation):

-25 °C to 60 °C (-13 °F to 140 °F)

Permissible ambient temperature (storage/transport):

-25 °C to 60 °C (-13 °F to 140 °F)

4 Product Specification

4.1 Simulator, Pre-registered Values

Order No.	Product Designation	Legend	pH Value	Temp.
SIM pH Kit			Order No. 52 300 410	
52 300 411	SIM PH4	SIM PH4	pH 4.00	80 °C
52 300 412	SIM PH7	SIM PH7	pH 7.00	25 °C
52 300 413	SIM PH TGL	Toggle Mode 30 sec	pH 4.00 & pH 7.00	25 °C & 80 °C
52 300 414	SIM PH ERR1	Glass Breakage	pH 7.00	25 °C
52 300 415	SIM PH ERR2	Fail Temperature Range	pH 10.67	-28 °C

4.2 Sensor Wear Data

Sensor Wear Monitor pH	Transmitters M400/M420/M700	% BAR
Time to maintenance	0 d	100 %
Life time indicator	245 d	18 %
Adaptive Cal Timer	7 d	0 %
Sensor operating time	55 of 300 d	
Autoclaving cycles	1 of 10	
CIP cycles	7 of 10	
SIP cycles	3 of 10	
Max temperature	130 °C dd.mm.yyyy (date of production)	

Exceptions

52 300 414 (SIM pH ERR1) DLI expired

52 300 415 (SIM pH ERR2) DLI expired

4.3 Warnings and Alarms

Warnings and Alarms**Transmitters M300/M400**

52 300 414 (SIM pH ERR1)	Glass Breakage
--------------------------	----------------

52 300 415 (SIM pH ERR2)	Fail Temperature Range
--------------------------	------------------------

Transmitter M420

52 300 414 (SIM pH ERR1)	Glass Breakage (Err 15)
--------------------------	-------------------------

52 300 415 (SIM pH ERR2)	Fail Temperature Range (Err 13)
--------------------------	---------------------------------

Transmitter M700

52 300 414 (SIM pH ERR1)	Glass Breakage
--------------------------	----------------


52 300 415 (SIM pH ERR2)	Fail Temperature Range
--------------------------	------------------------

Please refer to your transmitter instruction manual for specific instructions.

5 Installation

5.1 Dimensions

Markings

The markings labeled on each pH ISM plug-in head state the following information:

6 Ordering Information

Order No.	Designation
52 300 410	SIM pH Kit
52 300 411	SIM pH4
52 300 412	SIM pH7
52 300 413	SIM pH TGL
52 300 414	SIM pH ERR1
52 300 415	SIM pH ERR2
52 300 416	SIM O₂ Kit
52 300 417	SIM O ₂ ZERO
52 300 418	SIM O ₂ AIR
52 300 419	SIM O ₂ TGL
52 300 420	SIM O ₂ ERR1
52 300 421	SIM O ₂ ERR2
52 300 422	SIM O₂ ppb Kit
52 300 423	SIM O ₂ ppb ZERO
52 300 424	SIM O ₂ ppb AIR
52 300 425	SIM O ₂ ppb TGL
52 300 426	SIM O ₂ ppb ERR1
52 300 427	SIM O ₂ ppb ERR2
52 300 428	SIM O₂ Trace Kit
52 300 429	SIM O ₂ Trace ZERO
52 300 430	SIM O ₂ Trace AIR
52 300 431	SIM O ₂ Trace TGL
52 300 432	SIM O ₂ Trace ERR1
52 300 433	SIM O ₂ Trace ERR2

The single simulators are not available for sale.

7 Certificates

ATEX/FM certification pending.

8 Protocol

Protocol for pH ISM Qualification Kit

Facility Name _____

Instrument Manufacturer _____

Date of Purchase _____

Transmitter Model _____

Transmitter Serial Nr _____

Simulator P/N	Simulator S/N	Designation	Valid until	Parameter/s (check all that apply for this instrument)

Date/Time	Acceptable? ⇒ Yes	Acceptable? ⇒ No	Name/Signature

Protocol for pH ISM Qualification Kit © Mettler-Toledo AG, MarCom, Im Hackacker 15, CH-8902 Urdorf, 03/2009 1

A copy of this form in A4 format is delivered with the pH ISM Qualification Kit.

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Protocole pour set de qualification ISM pH © Mettler-Toledo AG, MarCom, Im Hackacker 15, CH-8902 Urdorf. 03/2009 1
- CH Mettler-Toledo (Schweiz) GmbH**, Im Langacher, Postfach, CH-8606 Greifensee
- CN Mettler-Toledo Instruments (Shanghai) Co. Ltd.**, 589 Gui Ping Road Cao He Jing,
- D Mettler-Toledo GmbH**, Prozeßanalytik, Ockerweg 3, D-35396 Gießen
- F Mettler-Toledo, Analyse Industrielle S.A.S.**, 30, Boulevard de Douaumont, BP 949,
- UK Mettler-Toledo LTD**, 64 Boston Road Beaumont Leys, GB-Leicester LE4 1AW