

# Excellence Weighing Solutions



## ErgoClips Unique Safe Dosing

ErgoClips allow easy weighing directly into tare containers



**Don't waste your time**  
Say "Goodbye" to weighing paper. Its a thing of the past.



**Perfect Positioning**  
ErgoClips guarantee secure placement of different tare containers.



**Safe Dosing**  
Dosing is fast and safe thanks to the adjustable MinWeigh Door.




**Accurate Weighing**  
With ErgoClips you can rely on complete sample transfer.

**METTLER TOLEDO**

## ErgoClips Overview





Micro  
Balances  
XP26  
and XP56




<b>Part</b>	ErgoClip Basket micro	ErgoClip Flask micro	MinWeigh Door micro	Stand for ErgoClips micro
<b>Size</b>	Ø max. 23 mm	1-10 ml	Inner draft shield with a small adjustable opening	Holder to store up to 3 ErgoClips
<b>Tare containers (examples)</b>	Vials, test tubes, centrifuge tubes	Volumetric flasks		
<b>Order No.</b>	11107889	11107879	11107869	11140175

## XS/XP Analytical Balances



<b>Part</b>	ErgoClip Weighing Boat	ErgoClip Titration Basket	ErgoClip Basket small	ErgoClip Tube	ErgoClip Filter Holder	ErgoClip Quantos and vial adapters
<b>Size</b>	LxW: 60x45 mm	Ø max. 50 mm	LxW: 32x32 mm	Ø 8 / 12 mm	Ø 110-150 mm	Ø 10.3-21 mm
<b>Tare containers (examples)</b>	Weighing boats	Titration beakers, test tubes	Test tubes, centrifuge tubes or small plastic tubes	PCR tubes, Eppendorf® tubes, reaction tubes	Filters	HPLC tubes, vials
<b>Order No.</b>	11106748	11106883	11106747	11106784	11140185	11141570



<b>Part</b>	ErgoClip Round Bottom Flask	ErgoClip Flask	ErgoClip Small Flask	MinWeigh Door	Stand for ErgoClips
<b>Size</b>	Ø min. 46 mm	20-100 ml	1-10 ml	Draft shield with a small adjustable opening	Holder to store up to 3 ErgoClips
<b>Tare containers (examples)</b>	Round bottom flasks, beakers	Volumetric flasks	Volumetric flasks		
<b>Order No.</b>	11106746	11106764	11140180	11106749	11140170

**Mettler-Toledo AG**  
Laboratory & Weighing Technologies  
CH-8606 Greifensee, Switzerland  
Phone +41 44 944 22 11  
Fax +41 44 944 30 60

Subject to technical changes  
© 03/2009 Mettler-Toledo AG  
Printed in Switzerland 11796297  
Global MarCom Greifensee

[www.mt.com/ergoclips](http://www.mt.com/ergoclips)

For more information